

Cannula Ø 2,5 mm tip


Technical data sheet

1. Administrative information about the company

Date of issue : 13/02/24

1.1	Name : HYGITECH SAS	
1.2	Complete address : 55 avenue Foch 75116 PARIS – FRANCE	Phone: +33 1 48 01 32 89 e-mail : info@hygitech.fr Web site : www.hygitech.fr
1.3	Contact information Materiovigilance : Said BOUTOUGHMAS	Phone : +33 1 45 91 31 82 e-mail : qualite@hygitech.fr

2. Information about the equipment

2.1	Trade name: Cannula Ø 2,5 mm tip	
2.2	Class of MD : Class IIa European Directive: MDR 2017/745 According Annex n° XI Part A Notified body: 1639 SGS BELGIUM Date of first launch on the UE market: 2012_08 Manufacturer : Hygitech	
2.3	Product description : The extra fine tip opening ensures precise aspiration without blocking the field of view. Specially designed for use with microscopes and loupes. Each box contains 2 autoclavable adapters. <div data-bbox="338 1252 1434 1798">  </div>	

2.4	Packaging : <table border="1" data-bbox="344 322 1484 416"> <thead> <tr> <th>Reference</th><th>OU Order unit</th><th>CARTON</th><th>MOQ Minimum order quantity</th></tr> </thead> <tbody> <tr> <td>HY-12031</td><td>Box of 20 units</td><td>Carton of 30 boxes of 20 units</td><td>1 box of 20 units</td></tr> </tbody> </table> <p>Labelling : Unknown</p>	Reference	OU Order unit	CARTON	MOQ Minimum order quantity	HY-12031	Box of 20 units	Carton of 30 boxes of 20 units	1 box of 20 units										
Reference	OU Order unit	CARTON	MOQ Minimum order quantity																
HY-12031	Box of 20 units	Carton of 30 boxes of 20 units	1 box of 20 units																
2.8	Composition of the product, dimensions and accessories : <table border="1" data-bbox="344 584 1353 669"> <thead> <tr> <th>Elements</th><th>Material</th><th>Dimensions</th></tr> </thead> <tbody> <tr> <td>Cannula</td><td>PVC</td><td>Ø 2,5 mm and L : 18 cm</td></tr> </tbody> </table> <table border="1" data-bbox="344 736 1484 916"> <thead> <tr> <th colspan="3">Presence of :</th></tr> </thead> <tbody> <tr> <td>Latex</td><td><input type="checkbox"/> YES</td><td><input checked="" type="checkbox"/> NO</td></tr> <tr> <td>Phthalates (DEHP)</td><td><input type="checkbox"/> YES</td><td><input checked="" type="checkbox"/> NO</td></tr> <tr> <td>Products of animal or biological origin</td><td><input type="checkbox"/> YES</td><td><input checked="" type="checkbox"/> NO</td></tr> </tbody> </table>	Elements	Material	Dimensions	Cannula	PVC	Ø 2,5 mm and L : 18 cm	Presence of :			Latex	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO	Phthalates (DEHP)	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO	Products of animal or biological origin	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO
Elements	Material	Dimensions																	
Cannula	PVC	Ø 2,5 mm and L : 18 cm																	
Presence of :																			
Latex	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																	
Phthalates (DEHP)	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																	
Products of animal or biological origin	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																	
2.9	Field - Indications: Unknown																		
3. Method of sterilizing																			
	Sterile : YES Sterilization method: Ethylene Oxide (EO)																		
4. Storage conditions																			
	Normal conditions of storage: Store in a dry environment and at room temperature. Special precautions: Keep this product away from direct sunlight. Shelf life: 3 years Presence of temperature indicators: NO																		
5. Safety in use																			
5.1	Technical security: N/A																		
5.2	Biosecurity: N/A																		
6. Advices of use																			
	Contents sterile - Unless this package is opened or damaged - Do not re-sterilise - Destroy after use																		
7. Additional product information																			
	Single use: YES																		
8. Related documents																			
	Unknown																		